

Tandem Affinity Purification (TAP) Protocol

Material: Rice suspension cultured cells

Protein Extract Preparation

- 1, Rice suspension cells were harvested in liquid nitrogen.
- 2, Frozen tissues were homogenized to a fine powder with a motor and pestle under evaporating liquid nitrogen.
- 3, The resultant powder was resuspended in 2 volumes of chilled TBSN (20mM sodium phosphate pH7.5, 300mM NaCl, 0.1% NP-40).
- 4, The mixture was filtered through 2 layers of Miracloth, centrifuged at 10,000g for 20min at 4°C, and filtered through a 0.45µm membrane (Millipore).

Purification and Elution of Protein

- 1, A 50mg portion of protein extract was diluted to 45ml in TBSN and incubated with 300µl IgG sepharose 6 fast flow (GE Healthcare) on a rotator at 4°C overnight.
- 2, After washing the beads 3 times with 10ml of TBSN, pour the beads to chromatography empty column (BioRad Poly Prep).
- 3, Wash the beads 3 times with 1ml of TEV cleavage buffer (50mM Tris-HCl, pH8.0, 0.5mM EDTA, 1mM DTT).
- 4, Close the bottom, and add 1ml of TEV cleavage buffer and 5µg of AcTEV protease (Invitrogen). The column was rotated at 4°C overnight.
- 5, Drain eluate into a new column sealed at the bottom.
- 6, Wash out old column with 1ml TEV cleavage buffer and mix eluate to the new column.
- 7, Add 3 volumn (6ml) of Calmodulin-binding buffer (25mM Tris-HCl, pH8.0, 150mM NaCl, 1mM Mg acetate, 1mM imidazole, 2mM CaCl₂, 10mM 2-Me) with 0.02% NP-40 and 6µl of 1M CaCl₂.
- 8, Add 300µl of calmodulin Sepharose 4B (GE Healthcare) and incubate on a rotary shaker for 1h at 4°C.
- 9, Wash twice with 1ml of Calmodulin-binding buffer with 0.1% NP-40.
- 10, Wash once with 1ml of Calmodulin-binding buffer with 0.02% NP-40.

11, Plug the bottom of the column and add 2ml of Calmodulin elution buffer (25mM Tris-HCl, pH8.0, 150mM NaCl, 0.02% NP-40, 1mM Mg acetate, 1mM imidazole, 20mM EGTA, 10mM 2-Me).