

Transient Plant Cell Transformation by Particle Bombardment

(from Wong Hann Ling Sept 2006)

<Reagents>

2.5 M CaCl₂ (Filtered, store at -20°C)

0.1 M Spermidine (Filtered, store at -20°C)

0.8% Agarose plates (Can be substituted with wet filterpaper in Petri dishes)

<Gold particle preparation>

1) Weigh gold particles

0.6 mg particle (diam. 1.6 µm for onion and rice, 1.0 µm for Arabidopsis)/1 shot

2) Sterilize gold particles

Add 300 µl 70% EtOH into weighed gold particles and vortex.

Briefly spin down and discard EtOH.

Add 300 µl of 100% EtOH, spin down and discard EtOH.

3) Coating DNA onto gold particles

(NOTE: Following is the most important step in this experiment. Prepare all the solution suck in a tip respectively, and add them one by one seamlessly)

Add sterile dH₂O and briefly sonicate.

(NOTE: total volume of water and DNA together is 10 µl)

Add 5 µg/shot of DNA, vortex briefly.

Add 10 µl/shot of 2.5 M CaCl₂ solution and vortex briefly.

Add 4 µl/shot of 0.1 M Spermidine and vortex briefly.

Incubate at room temp. for 30 min.

Spin down and discard supernatant.

Add 100 µl 70%EtOH and sonicate briefly, spin down and remove supernatant.

Add 15 µl 100%EtOH/ shot and briefly sonicate.

Spot resuspended samples on microcarrier and allow to them dry.

<Operation particle gun (BioRad PDS-1000/He)>

- 1) Turn on Vacuum pump.
- 2) Turn on He meter gauge.
- 3) Turn on He meter valve.
- 4) Turn the gauge to the left of He meter and set it at 1300 psi (depending on plant tissue and rapture disk types)
- 5) Place rapture disk (e.g. 1100 psi type) on Gas Acceleration Tube Stand and tighten with Torque wrench.
- 6) Place Stopping Screen and Microcarrier on Microcarrier Launch Assembly in main chamber (at 2nd level).
- 7) Place plant tissue on agarose plate and set them on the 4th level in main chamber.
- 8) Shut the door of chamber apply vacuum by pressing <VAC> down until vacuum reached 28.5~29.0, then quickly press down <HOLD>.
- 9) Press down <FIRE> to launch bombardment. (A 'ponk' sound with be emitted when bombardment occurred.)
- 10) Release vacuum by placing switching knob at <VENT> position.
- 11) Repeat steps 4~8 if required.
- 12) To finish, turn off he valve, shut the cahmber, and set <VAC> and <FIRE> simultaneously until vacuum He pressure returned zero.
- 13) Turn off He meter valve.
- 14) Add a little water to bombarded plates and store in dark at 30°C.

*2 DNA amount is depending on the plant. If you are using Arabidopsis or onion, you may use 0.5~1 µg, but bombarding into rice requires more DNA